

	Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT	KAMSIDAR		Page 1 of 72
BATCH SIZE	600,000 Tablets	Mfg. Date :	
BATCH No.		Exp. Date :	
Effective Date:	MFR No.: KPI/MFR/022/00		

Generic Name:	Sulfadoxine and Pyrimethamine Tablets USP 525mg		BMR Revision No: 03	Approvals	Start date
Product Pack No:(√/ mark as appropriate)	KSD1TB	KSD2TJ	KSD3TJ	QA Head	
Pack Size: (√/ mark as appropriate)	10x10 Blisters	500's	1000's	Company/Prod. Pharmacist	
Reference SOP Number:	QAD/022/01			Production Manager	

RAW MATERIAL REQUISITION/DISPENSING SHEETS

Sr.No.	Description	A.R. No.***	Unit Quantity (mg/tablet)	Overage (mg/tablet)	Total Quantity (mg/tablet) *	Total Batch Quantity (Kg)	Lot	Quantity Per Lot	Tare Weight	Gross Weight	Net Weight	Scale ID	Weighed by	Checked by**	Received by
1	Sulfadoxine BP SDX1XN		500.000	0.000	500.000	300.000	1	75.00kg							
							2	75.00kg							
							3	75.00kg							
							4	75.00kg							
2	Microcrystalline Cellulose BP MCC2XN		40.000	0.000	40.000	24.000	1	6.00kg							
							2	6.00kg							
							3	6.00kg							
							4	6.00kg							

PREPARED BY: Sign & Date:	REVIEWED BY: Sign & Date:	APPROVED BY: Sign & Date:	AUTHORIZED BY : Sign & Date:
Name:	Name:	Name:	Name:
Designation:	Designation:	Designation:	Designation:

	Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT	KAMSIDAR		Page 2 of 72
BATCH SIZE	600,000 Tablets	Mfg. Date :	
BATCH No.		Exp. Date :	
Effective Date:	MFR No.: KPI/MFR/022/00		

Generic Name:	Sulfadoxine and Pyrimethamine Tablets USP 525mg		BMR Revision No: 03	Approvals	Start date
Product Pack No:(√/ mark as appropriate)	KSD1TB	KSD2TJ	KSD3TJ	QA Head	
Pack Size: (√/ mark as appropriate)	10x10 Blisters	500's	1000's	Company/Prod. Pharmacist	
Reference SOP Number:	QAD/022/01			Production Manager	

RAW MATERIAL REQUISITION/DISPENSING SHEETS

Sr.No.	Description	A.R. No.***	Unit Quantity (mg/tablet)	Overage (mg/tablet)	Total Quantity (mg/tablet) *	Total Batch Quantity (Kg)	Lot	Quantity Per Lot	Tare Weight	Gross Weight	Net Weight	Scale ID	Weighed by	Checked by**	Received by
3	Maize Starch BP MZS2XN (dry mix)		40.000	0.000	40.000	24.000	1	6.00kg							
							2	6.00kg							
							3	6.00kg							
							4	6.00kg							
4	Maize Starch BP MZS2XN (For paste)		20.000	0.000	20.000	12.000	1	3.00kg							
							2	3.00kg							
							3	3.00kg							
							4	3.00kg							

PREPARED BY: Sign & Date:	REVIEWED BY: Sign & Date:	APPROVED BY: Sign & Date:	AUTHORIZED BY : Sign & Date:
Name:	Name:	Name:	Name:
Designation:	Designation:	Designation:	Designation:

	Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT	KAMSIDAR		Page 3 of 72
BATCH SIZE	600,000 Tablets	Mfg. Date :	
BATCH No.		Exp. Date :	
Effective Date:	MFR No.: KPI/MFR/022/00		

Generic Name:	Sulfadoxine and Pyrimethamine Tablets USP 525mg		BMR Revision No: 03	Approvals	Start date
Product Pack No:(√/ mark as appropriate)	KSD1TB	KSD2TJ	KSD3TJ	QA Head	
Pack Size: (√/ mark as appropriate)	10x10 Blisters	500's	1000's	Company/Prod. Pharmacist	
Reference SOP Number:	QAD/022/01			Production Manager	

RAW MATERIAL REQUISITION/DISPENSING SHEETS

Sr.No.	Description	A.R. No.***	Unit Quantity (mg/tablet)	Overage (mg/tablet)	Total Quantity (mg/tablet) *	Total Batch Quantity (Kg)	Lot	Quantity Per Lot	Tare Weight	Gross Weight	Net Weight	Scale ID	Weighed by	Checked by**	Received by
5	Propylparaben BP PPB2XN		0.180	0.000	0.180	0.108	1	27.00g							
							2	27.00g							
							3	27.00g							
							4	27.00g							
6	Methylparaben BP MPB2XN		1.230	0.000	1.230	0.738	1	184.5g							
							2	184.5g							
							3	184.5g							
							4	184.5g							
7	Pyrimethamine BP PYM1XN		25.00	0.00	25.00	15.00	1	15.00kg							

PREPARED BY: Sign & Date:	REVIEWED BY: Sign & Date:	APPROVED BY: Sign & Date:	AUTHORIZED BY : Sign & Date:
Name:	Name:	Name:	Name:
Designation:	Designation:	Designation:	Designation:

	Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT	KAMSIDAR		Page 4 of 72
BATCH SIZE	600,000 Tablets	Mfg. Date :	
BATCH No.		Exp. Date :	
Effective Date:	MFR No.: KPI/MFR/022/00		

Generic Name:	Sulfadoxine and Pyrimethamine Tablets USP 525mg		BMR Revision No: 03	Approvals	Start date
Product Pack No:(√/ mark as appropriate)	KSD1TB	KSD2TJ	KSD3TJ	QA Head	
Pack Size: (√/ mark as appropriate)	10x10 Blisters	500's	1000's	Company/Prod. Pharmacist	
Reference SOP Number:	QAD/022/01			Production Manager	

RAW MATERIAL REQUISITION/DISPENSING SHEETS

Sr.No.	Description	A.R. No.***	Unit Quantity (mg/tablet)	Overage (mg/tablet)	Total Quantity (mg/tablet) *	Total Batch Quantity (Kg)	Lot	Quantity Per Lot	Tare Weight	Gross Weight	Net Weight	Scale ID	Weighed by	Checked by**	Received by
8	Magnesium Stearate BP MAG2SN		6.00	0.00	6.00	3.60	1	3.60kg							
9	Sodium Lauryl Sulphate BP SLS2XN		8.31	0.00	8.31	4.99	1	4.99kg							
10	Sodium Starch Glycollate BP SSG2XN		20.00	0.00	20.00	12.00	1	12.00kg							
11	Colloidal Anhydrous Silica BP (Aerosil 200) AER2XN		13.86	0.00	13.86	8.32	1	8.32kg							
TOTAL			674.58	0.00	674.58	404.75				Store-in charge	Dispensing SPV		QA Officer		
Weighing Started (HH:MM)			Date:		Weighing Stopped (HH:MM)			Date:							
										(Sign & Date)	(Sign & Date)		(Sign & Date)		

*Includes conversion to base drug, with overage and corrected for Potency and/or Moisture where required.

**Check weight, A.R. No. and verify presence of QC approved label on Raw Material container.

***AR No. Checked/Approved by QCM _____
(Sign & Date)

PREPARED BY: Sign & Date:	REVIEWED BY: Sign & Date:	APPROVED BY: Sign & Date:	AUTHORIZED BY : Sign & Date:
Name:	Name:	Name:	Name:
Designation:	Designation:	Designation:	Designation:

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 5 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
PRODUCT DETAILS				

Label Claim	:	Each Kamsidar Tablet Contains: Sulfadoxine and Pyrimethamine Tablets USP 525mg Colour: White
Shelf Life	:	3 years
Storage Conditions	:	Store in a cool & dry place. Protect from direct light, heat and moisture.
Mfg. License. Nos.	:	NDA/MAL/HDP/5467
Brand Name /Product Name	:	Kamsidar
Required Batch Size	:	600,000 Tablets
Market	:	Domestic/Local
Average weight (Uncoated)	:	674.58mg
Dissolution	:	Not Less than 65 % for each Tablet
Assay:	:	Percentage Limit (90.0 % to 110.0 % of the Label claim)

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 6 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
INDEX				

Sr. Nos.	TITLES	Page Nos.
1.	Raw Material Requisition/Dispensing Sheets	1-3
2.	Product Details	4
3.	Index	5
4.	Abbreviations	6
5.	Equipment List& Instruments List	7
6.	General Instruction	8
7.	Safety Precaution	8
8.	General SOP In Manufacturing Respectively	8
9.	Dispensing Line Clearance	9
10.	Mixing/Granulation	10-15
11.	Blending/Lubrication	16-19
12.	Slugging	20-25
13.	Re-Blending	26-29
14.	Tablet Compression	30-47
15.	Tablet Inspection & Sorting	48-49
16.	QA In-Process Quality Control Report For Inspection & Sorting	50
17.	Packaging	52-57
18.	In-Process Quality Control Report For Primary Packing	58-60
19.	In-Process Quality Control Report For Secondary Packing	63-65
20.	Finished Product Reconciliation Sheet	66
21.	Batch Release Sheet	67
22.	Revision History	68

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 7 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
ABBREVIATIONS				

S/N	SHORT	IN FULL	S/N	SHORT	IN FULL
01	BMR	Batch Manufacturing Record	19	L.O.D.	Loss On Drying
02	BPR	Batch Packing Record	20	NMT	Not More Than
03	API	Active Pharmaceutical Ingredients	21	NLT	Not Less Than
04	LAF	Laminar Air Flow	22	QAO	Quality Assurance Officer
05	RLAF	Reverse Laminar Air Flow	23	SPV	Supervisor
06	SOP	Standard Operating Procedure	24	IPQA	In-Process Quality Assurance
07	B/L	Break Line	25	IPA	Isopropyl Alcohol
08	AR No.	Analytical Report Number	26	NA	Not Applicable
09	Mfg.	Manufacturing	27	QA	Quality Assurance
10	Exp.	Expiry	28	M/C	Machine
11	RMG	Rapid Mixer Granulator	29	RPM	Revolutions Per Minute
12	FBD	Fluidized Bed Dryer	30	PMF No.	Product Master Formula Number
13	B. size	Batch Size	31	PVC	Polyvinyl Chloride
14	HDPE	High Density Polyethylene	32	BOPP	Biaxially Oriented Polypropylene
15	SMFPQ	Semi-Finished Product Quarantine	33	QC	Quality Control
16	IPQC	In-Process Quality Control	34	S/N	Serial Number
17	FIFO	First In-First Out	35	FEFO	First Expiry-First Out
18	LHS	Left Hand Side	36	RHS	Right Hand Side

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 8 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
LIST OF EQUIPMENTS & INSTRUMENTS				

List of Equipment used in the Manufacturing of Product.

S/N	Names	Mark	S/N	Names	Mark
DISPENSING					
1	Weighing Balance B-46		IPQC LAB		
2	Weighing Balance G-28				
3	Weighing Balance B-59				
4	RLAF Booth LAF-01		1	Weighing Balance B-53	
5	RLAF Booth LAF-05		2	Weighing Balance B-54	
GRANULATION 1			3	Weighing Balance B-55	
1	Rapid Mixer Granulator T-47		4	Weighing Balance B-23	
2	Fluidized Bed Dryers T-48		5	Weighing Balance B-52	
3	Fluidized Bed Dryers T-67		6	IR Moisture Analyzer B-19	
4	Vibro Sifter T-77		7	IR Moisture Analyzer B-43	
5	Paste Kettle T-49		8	Friability Apparatus QD-075	
6	Comminuting Mill T-50		9	Friability Apparatus T-33	
7	Comminuting Mill T-14		10	Hardness Tester QD-081	
8	Weighing Balance B-60		11	Hardness Tester QD-073	
GRANULATION 2			12	Disintegration Apparatus QD-084	
1	Rapid Mixer Granulator T-27		13	Disintegration Apparatus QD-085	
2	Fluidized Bed Dryers T-24		14	Vernier Caliper QD-086	
3	Fluidized Bed Dryers T-26		15	Vernier Caliper QD-087	
4	Vibro Sifter T-78		16	Leak Test Apparatus QD-088	
			17	Friability Apparatus QD-077	
5	Paste Kettle T-60		VISUAL INSPECTION AND SORTING		
6	Multimill T-44		1	Metal Detector	
7	Comminuting Mill T-01		2	De-duster	
8	Weighing Balance B-13		BLISTER PACKING		
COMPRESSION			1	Blister Machine P-03	
1	Compression Machine T-06		2	Blister Machine P-10	
2	Compression Machine T-15		3	Blister Machine P-20	
3	Compression Machine T-30		4	Blister Machine P-21	
4	Compression Machine T-35		5	Blister Machine P-27	
5	Compression Machine T-52		STRIP PACKING		
6	Compression Machine T-80		1	Strip Packing Machine	
BLENDING			OTHERS (SPECIFY)		
1	Double Cone Blender T-09				
2	Double Cone Blender T-17				
3	Double Cone Blender T-04				
4	Vibrosifter T-66				
5	Vibrosifter T-79				

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 9 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	

GENERAL INSTRUCTIONS

1. Use approved Equipment only, mentioned in respective stages of the BMR.
2. All Equipment for use must be duly cleaned as per respective SOPs and approved by QA.
3. Each container, machinery or Equipment must bear Status indicating label, duly signed by the person responsible.
4. All Equipment used, must have appropriate safety measures.
5. Care should be taken to ensure that all manufacturing Operations are performed as per the process under supervision of competent technical staff, workers should wear nose masks, and clean gloves, gowning throughout the manufacturing and packing operations.
6. Ensure that all the ingredients have Release labels for intact containers and proper A. R. Nos. for loose containers.
7. Ensure the materials are dispensed and checked as per respective SOPs and FIFO / FEFO system is followed.
8. Dispensing materials should be performed in presence of three qualified people, 1 from stores, 1 from Q.A. and 1 from production.
9. Observe the atmospheric conditions mentioned at various steps, in the BMR.
10. Dispensing labels shall be affixed with BMR.

SAFETY PRECAUTIONS

1. Safety Precautions to be taken to avoid explosions in FBD during drying.
2. Safety Precautions must be taken during grinding, milling, granulating process.
3. Protect Respiratory Organs by using cloth masks and gloves during manufacturing process.
4. Store all granules and tablets in well closed container lined with double polybags.
5. Follow personal hygienic requirements.
6. To avoid any accident, check the earthing and continuity of electricity while operating the machine.
7. Check the safety mechanisms at the start of the operations, at the end of operations, after power failures and any breakdowns of the machine.
8. **Material Transfers:** -Transfer all the dispensed materials from stores to manufacturing areas. Counter-check and record identities and weights of all the ingredients received. Do not use any material which does not comply with weight or analytical report (AR) numbers.

Storage of Materials: -Store the material, in process material, semi-finished and finished material in designated area properly labeled mentioning status of the material and batch details.

NOTE: Ensure all the SOPs used for references in the BMR are of current versions.

GENERAL SOPs IN MANUFACTURING RESPECTIVELY.		
S. Nos.	Titles	SOP Nos.
1.	Procedure for Environmental Monitoring Temperature & Differential Pressure.	QAD/040/00
2.	Procedure for Area Cleaning and sanitization of Premises.	GEN/008/06
3.	In-Process Quality Control and Assurance	QAD/021/00

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 10 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
DISPENSING					

LINE CLEARANCE FORM				
Previous product:		Batch No:		
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by
				QA (Sign & Date)
1. Area cleaning is done as per SOP.				
2. Ensure the absence of batch documents, labels, materials or remnants of previous product or batch.				
3. Ensure availability of the BMR for the in-coming batch and Verify that ALL materials and facilities for new batch are available, labeled and identified.				
4. Gowning and de-gowning procedure is followed.				
5. Dispensing Room, Dispensing booth, Balances and dispensing tools are clean and have status as ready for use.				
6. Clean scoops, containers and poly bags are available for dispensing as per SOP.				
7. Manometer reading of dispensing booth RLAF is within the required limit. Actual Manometer reading:				
8. Balance Calibration carried out before starting dispensing.				
9. Verify the release status & retest validity period of each material.				
10. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (40-65)% _____				
11. Ensure that the RLAF is started 15 minutes before start of the dispensing activities.				
Remarks:				
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by
				QA (Sign & Date)
1. All material from current product have been removed.				
2. All facilities have been cleaned and labeled.				
3. The area has been cleaned and labeled.				
4. All paperwork for the current batch has been completed.				
5. Reconciliation of current batch has been done.				
6. Ensure that utensils and accessories from previous operations have been removed.				
Remarks:				

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 11 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MIXING/GRANULATION					

LINE CLEARANCE FORM				
Previous product:		Batch No:		
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)
				Status
1. Area cleaning as per SOP.				
2. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.				
3. Gowning procedure is followed.				
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.				
5. Check the proper cleanliness of Vibro-sifter, Paste Kettle, RMG, FBD, Multi mill etc.				
6. Ensure that the Vibrosifter sieves, Multimill screen, FBD bags, Scoops etc. cleaned.				
7. Check and ensure that Balance Calibration & Verification records are updated.				
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (40-65)% _____				
9. Check the Integrity/Intactness of FBD gaskets, Product Bowl sieves and Dutch Meshes/Distribution Plates, and the sifting and milling sieves.				
Remarks:				
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)
				Status
1. All material from current product have been removed.				
2. Check the cleanliness of Supply and Return air Grills, Light Fixtures, Waste Bins, and drainage points.				
3. The area has been cleaned and labeled.				
4. All paperwork for the current batch has been completed.				
5. Reconciliation of current batch has been done.				
6. Machine(s) has been checked for cleanliness.				
Remarks:				

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 12 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
MIXING/GRANULATION				

There are four lots of material for this batch, for each perform the following:

A. DRY MIXING

- Pulverize each lot of Sulfadoxine and Pyrimethamine by passing through Multimill fitted with 0.5mm S.S. sieve, impact forward at speed 2.
- In clean stainless steel Rapid Mixer Granulator, add in sequence:
 - 75.00kg of Sulfadoxine**
 - 6.00kg of Maize Starch**
 - 6.00kg of Microcrystalline Cellulose**
- Dry blend at SLOW speed in Rapid Mixer Granulator for 10 minutes then at FAST speed for 5 minutes.

Lot	RMG NUMBER:									
	Speed	Impeller Speed [RPM]	Chopper Speed [RPM]	Started	Completed	Actual Time Taken (min)	Time Specified (min)	Done by	Verified by SPV	Confirmed by QA
1	Slow									
	Fast									
2	Slow									
	Fast									
3	Slow									
	Fast									
4	Slow									
	Fast									

B. WET MIXING

Prepare starch paste as follows:

- Boil 24.0 litres of purified Water in stainless steel vessel and add while mixing.
 - 27.00g of Propylparaben**
 - 184.5g of Methylparaben**
- Meanwhile, add in sequence in a clean stainless steel mixing vessel.
 - 6litres of cold purified water**
 - 3.00kg of Maize Starch**
- Mix until smooth uniform slurry forms.
- Add the slurry to the hot solution in 4 above while vigorously stirring to get a smooth uniform paste.
- Add the paste in 6 to the dry blend in 3 and mix at SLOW speed for 5 minutes then at FAST speed for 25 minutes.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 13 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MIXING/GRANULATION					

Lot	RMG NUMBER:														
	Binder Preparation				Massing/Mixing										
	STD Cold Water Volume	Actual Cold Water Volume	STD Hot Water Volume	Actual Hot Water Volume	Speed	Impeller Speed [RPM]	Chopper Speed [RPM]	Started	Completed	Actual Time Taken (min)	Time Specified (min)	Done by	SPV	QA	
1					Slow										
					Fast										
2					Slow										
					Fast										
3					Slow										
					Fast										
4					Slow										
					Fast										

- Dry the granules from 8 above at 70 °C for 45 minutes in Fluid Bed Drier then remove and turn around the granules with a stainless steel scoop. Repeat the above procedure then mill the granules in Comminuting Mill fitted with sieve size 4. Redry the granules at 70 °C for 45 minutes.

9. Repeat the above procedure for the rest of the lots.

FIRST DRYING

LOT	FBD No.	DRYING CONDITIONS			DRYING TIME				Done by	SPV	QA
		STD Inlet Temp. (°C)	Observed Drying Temp. (°C)	Outlet Temp. (°C)	Started	Completed	Actual Time Taken (min)	Time Specified (min)			
1											
2											
3											
4											

SECOND DRYING

LOT	FBD No.	DRYING CONDITIONS			DRYING TIME				Done by	SPV	QA
		STD Inlet Temp. (°C)	Observed Drying Temp. (°C)	Outlet Temp. (°C)	Started	Completed	Actual Time Taken (min)	Time Specified (min)			
1											
2											
3											
4											

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 14 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MIXING/GRANULATION					

FINAL DRYING														
Specifications: LOD Limits: _____ %.														
L O T	FBD NO.	DRYING CONDITIONS			DRYING TIME				FINAL MILLING					
		STD Inlet Temp.(°C)	Observed Drying Temp.(°C)	Outlet Temp.(°C)	Started	Completed	Actual Time Taken (min)	Time Specified (min)	Miller No.	Sieve Size [mm]	LOD (%)	Done by	SPV	QA
1														
2														
3														
4														

10. Let the granules cool to room temperature then sift through #14 on vibratory sifter then mill granules with Multimill fitted with 3mm sieve with knives forward at speed 2.

QA REPORT

11. Loss on Drying: Check the moisture content in the granules:

Sample: Take 20g of dried granules and operate the moisture balance as per SOP.

Specifications: LOD Limits: _____ %.

Acceptance criteria: % LOD should comply with specification:

Test Time	Lot no.	Drying Temp.(°C)	Drying duration (min)	FBD no.	LOD (%)	Moisture Balance ID	Remarks	Done by QA (Sign & Date)

Attach the slip from the moisture analyzer.

Release procedure: Yes / No

If the above tests fall within the specified limits, QA fills out the results in the above table and the material is sent to the next processing stage.

Rejection procedure: Yes / No

If the material fails any of the above tests, QA fills out and labels the RMG bowl with "REJECTED / REWORK" label and the material is sent for reworking.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 15 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MIXING/GRANULATION					

12. Determine the percentage yield of the dried granules.

Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
Total			
Theoretical Weight			

Operator/ grp. Leader (sign & date) _____ Section Supervisor (sign & date) _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 16 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MIXING/GRANULATION					

YIELD RECONCILIATION				
Steps	Descriptions	Qty. (Kg)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size with respect to dispensed material			
B.	Actual Quantities of Dried Granules			
C.	Sample Quantities (L.O.D/In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
Percentage yield = $F \times 100 / A = \dots\dots\dots\%$				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):
CAUSE OF VARIATION: YES/NO.

REMARKS/BREAKDOWNS.....
.....
.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 17 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
BLENDING/LUBRICATION					

LINE CLEARANCE FORM					
Previous product:		Batch No:			
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. Area cleaning as per SOP.					
2. Check out and ensure that all the previous product containers, materials, and labels are removed from the manufacturing area.					
3. Gowning procedure is followed.					
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.					
5. Check the proper cleanliness of the Air Supply and Return grills.					
6. Ensure that the Vibrosifter sieves, Double Cone Blender, Scoops etc. are cleaned.					
7. Check and ensure that Balance Calibration & Verification records are updated.					
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (40-65)% _____					
Remarks:					
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. All material from current product have been removed.					
2. All facilities have been cleaned and labeled.					
3. The area has been cleaned and labeled.					
4. All paperwork for the current batch has been completed.					
5. Reconciliation of current batch has been done.					
6. Machine(s) has been checked for cleanliness.					
Remarks:					

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 18 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
BLENDING/LUBRICATION					

1. Equipment used for Sifting and Blending

Equipment Name	Equipment ID No.	Mark (Equipment used)
Double Cone Blender	T-04	
Double Cone Blender	T-09	
Double Cone Blender	T-17	
Vibrosifter	T-66	
Vibrosifter	T-79	

2. Sift using Vibrosifter with sieves indicated in the table below then transfer into Double Cone Blender

- _____ **Kg of Dried Granules.**
- 12.00kg Sodium Starch Glycollate**
- 3.60kg Magnesium Stearate**
- 4.99kg Sodium Lauryl Sulfate**
- 15.00kg Pyrimethamine**
- 8.32kg Aerosil 200**

3. Add the Recoveries, if any (indicate quantities and batch Numbers)

Batch No.				
Quantity				

S/N	Ingredient	Sift through mesh size	Qty added (Kg)	Added by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)
1	Dried Granules					
2	Sodium Starch Glycollate	#40				
3	Magnesium Stearate	#40				
4	Sodium Lauryl Sulfate	#40				
5	Pyrimethamine	#40				
6	Aerosil 200	#40				
8	Recoveries					
Total Weight of dry ingredients added						

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 19 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
BLENDING/LUBRICATION					

4. After all ingredients have been added to the double cone blender, mix for 15 minutes.

Start time	Completion time	Time taken	specifications	Deviation	Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)

NB: If recoveries are added, Calculate the new batch size:

Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)

5. Determine the percentage yield of the blended granules.

Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
Total			
Theoretical Weight			

Operator/ grp. Leader (sign & date) _____ Section Supervisor (sign & date) _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 20 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
BLENDING/LUBRICATION					

YIELD RECONCILIATION				
Steps	Descriptions	Qty. (Kg)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size with respect to Dried Granules + Blending Materials			
B.	Actual Quantities of Lubricated Granules			
C.	Sample Quantities (L.O.D/In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
Percentage yield = $F \times 100 / A = \dots\dots\dots\%$				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):
CAUSE OF VARIATION:

REMARKS/BREAKDOWNS.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 21 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SLUGGING					

LINE CLEARANCE FORM					
Previous product:			Batch No:		
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)	Status
1. Area cleaning as per SOP.					
2. Check and ensure that all the previous product containers, material, and labels are removed from the manufacturing area.					
3. Gowning procedure is followed.					
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp. Date & Status with sign & date.					
5. Ensure the absence of batch documents, labels, materials or remnants of previous product or batch on each critical part.					
6. Are the following parts of the compression machine cleaned: Hopper/Feeder/Turret/Below the turret, Y-chute/Powder collection, Punch & Die Set.					
7. Check and ensure that Balance Calibration & Verification records are updated.					
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (55 ± 10%) _____					
Remarks:					
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)	Status
1. All material from current product have been removed.					
2. All facilities have been cleaned and labeled.					
3. The area has been cleaned and labeled.					
4. All paperwork for the current batch has been completed.					
5. Reconciliation of current batch has been done.					
6. Machine(s) has been checked.					
Remarks:					

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 22 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SLUGGING					

IN-PROCESS QUALITY CONTROL CHECKS DURING SLUGGING									
Done by:		Operator:		QA:		Operator:		QA:	
Time									
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance	White round slugs without any colour spots (Black or Brown)								
Hardness	4-8 (kg/cm ²)								

IN-PROCESS QUALITY CONTROL CHECKS DURING SLUGGING									
Done by:		Operator:		QA:		Operator:		QA:	
Time									
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance	White round slugs without any colour spots (Black or Brown)								
Hardness	4-8 (kg/cm ²)								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 23 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SLUGGING					

IN-PROCESS QUALITY CONTROL CHECKS DURING SLUGGING									
Done by:		Operator:		QA:		Operator:		QA:	
Time									
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance	White round slugs without any colour spots (Black or Brown)								
Hardness	4-8 (kg/cm ²)								

IN-PROCESS QUALITY CONTROL CHECKS DURING SLUGGING									
Done by:		Operator:		QA:		Operator:		QA:	
Time									
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance	White round slugs without any colour spots (Black or Brown)								
Hardness	4-8 (kg/cm ²)								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 24 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SLUGGING					

Make sure the line clearance is carried out as per SOP for line clearance and recorded in the logbook before Slugging process.

1. Take the blended material for slugging.

Started	Date:	Time:	Completed	Date:	Time:
---------	-------	-------	-----------	-------	-------

2. Determine the percentage yield of the Slugs.

Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
Total			
Theoretical Weight			

Operator/ grp. Leader (sign & date) _____ Section Supervisor (sign & date) _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 25 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SLUGGING					

YIELD RECONCILIATION				
Steps	Descriptions	Qty. (Kg)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size with respect to Dried Granules + Blending Materials			
B.	Actual Quantities of Slugs			
C.	Sample Quantities (L.O.D/In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
Percentage yield =F x100/A =.....%				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):
CAUSE OF VARIATION:

REMARKS/BREAKDOWNS.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 26 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MILLING THE SLUGS					

LINE CLEARANCE FORM					
Previous product:		Batch No:			
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. Area cleaning as per SOP.					
2. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.					
3. Gowning procedure is followed.					
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp. Date & Status with sign & date.					
5. Ensure the absence of batch documents, labels, materials or remnants of previous product or batch on each critical part.					
6. Are the following parts of the Multimill machine cleaned: Milling screens/Sieves, the swing beaters/Blades.					
7. Check and ensure that Balance Calibration & Verification records are updated.					
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (55 ± 10%) _____					
Remarks:					
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. All material from current product have been removed.					
7. All facilities have been cleaned and labeled.					
8. The area has been cleaned and labeled.					
9. All paperwork for the current batch has been completed.					
10. Reconciliation of current batch has been done.					
11. Machine(s) has been checked for cleanliness.					
Remarks:					

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 27 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MILLING THE SLUGS					

1. After slugging, mill the material with Multimill No. _____, fitted with sieve # 3 and #4.

Started	Date:	Time:	Completed	Date:	Time:
---------	-------	-------	-----------	-------	-------

2. Determine the percentage yield of the Re-sized granules.

Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
Total			
Theoretical Weight			

Operator/ grp. Leader (sign & date) _____ Section Supervisor (sign & date) _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 28 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
MILLING THE SLUGS					

YIELD RECONCILIATION				
Steps	Descriptions	Qty. (Kg)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size with respect to slugs			
B.	Actual Quantities of Re-sized granules			
C.	Sample Quantities (L.O.D/In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
$\% \text{ yield} = F \times 100 / A = \dots\dots\dots\%$				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):
CAUSE OF VARIATION:

REMARKS/BREAKDOWNS.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 29 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
RE-BLENDING					

LINE CLEARANCE FORM				
Previous product:			Batch No:	
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.			Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)
			Verified by QA (Sign & Date)	Status
1. Area cleaning as per SOP.				
2. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.				
3. Gowning procedure is followed.				
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.				
5. Check the proper cleanliness of the Air Supply and Return grills.				
6. Ensure that the Vibrosifter sieves, Double Cone Blender, Scoops etc. are cleaned.				
7. Check and ensure that Balance Calibration & Verification records are updated.				
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (55 ± 10%) _____				
Remarks:				
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements			Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)
			Verified by QA (Sign & Date)	Status
1. All material from current product have been removed.				
2. All facilities have been cleaned and labeled.				
3. The area has been cleaned and labeled.				
4. All paperwork for the current batch has been completed.				
5. Reconciliation of current batch has been done.				
6. Machine(s) has been checked.				
Remarks:				

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 30 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
RE-BLENDING					

1. Re-blend the material in the double cone blender No. _____, mix for 20 minutes.

Start time	Completion time	Time taken	specifications	Deviation	Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date)

Intermediate Sample for moisture content and homogeneity of blended / lubricated granules.

1. Take samples of blended granules as per SOP for sampling blended / lubricated granules and send them to QA lab along with this record.

Quantity sampled	Top.....g	Middle.....g	Bottom.....g
Sampled by QA			

QA REPORT				
Test	Specification	Top Sample	Middle Sample	Bottom Sample
LOD by IR moisture balance (%)	1.0-3.0			

The samples **Comply** ☐ **Do not comply** ☐ With the above specifications.

QA Officer: (sign)..... Date.....

Release procedure: Yes / No

If the above tests fall within the specified limits, QA fills out and labels the blender with "PASSED" label and the material is sent to the next processing stage.

Rejection procedure: Yes / No

If the material fails any of the above tests, QA fills out and labels the blender with "REJECTED / REWORK" label and the material is sent for reworking.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 31 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
RE-BLENDING					

2. Determine the percentage yield of the blended granules.

Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
Total			
Theoretical Weight			

Operator/ grp. Leader (sign & date) _____ Section Supervisor (sign & date) _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 32 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
RE-BLENDING					

YIELD RECONCILIATION				
Steps	Descriptions	Qty. (Kg)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size with respect to Re-sized granules			
B.	Actual Quantities of blended Granules			
C.	Sample Quantities (L.O.D/In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
% yield =F x100/A =.....%				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):
CAUSE OF VARIATION:

REMARKS/BREAKDOWNS.....
.....
.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 33 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

LINE CLEARANCE FORM					
Previous product:		Batch No:			
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. Area cleaning as per SOP.					
2. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.					
3. Gowning procedure is followed.					
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.					
5. Ensure the absence of batch documents, labels, materials or remnants of previous product or batch on each critical part.					
6. Are the following parts of the compression machine cleaned: Hopper/Feeder/Turret/Below the turret, Y-chute/Powder collection, Punch & Die Set.					
7. Check and ensure that Balance Calibration & Verification records are updated.					
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (55 ± 10%) _____					
Remarks:					
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. All material from current product have been removed.					
2. All facilities have been cleaned and labeled.					
3. The area has been cleaned and labeled.					
4. All paperwork for the current batch has been completed.					
5. Reconciliation of current batch has been done.					
6. Machine(s) has been checked for cleanliness.					
Remarks:					

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 34 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

1. EQUIPMENT USED FOR COMPRESSION

Equipment Name	Equipment ID	Mark	Equipment Name	Equipment ID	Mark
Compression Machine	T-06		Compression Machine	T-80	
Compression Machine	T-15		Metal Detector	T-84	
Compression Machine	T-30		Metal Detector	T-85	
Compression Machine	T-35		Deduster	T-86	
Compression Machine	T-52		Deduster	T-87	

2. MACHINE SET UP

TYPE OF TOOLING :-D-TOOLING		B-TOOLING	
Punch size: (for Trade)	12.5mm	Punch shape: (for Trade).	
Tablet face appearance/ description (Trade): White, round, flat, beveled tablets embossed " F" on one side and double cross line on the other side.		Upper:	Flat, beveled, embossed "F"
		Lower:	Flat, beveled with double cross line
Punch size: (for UG)	12.5mm	Punch shape: (for UG).	
Tablet face appearance/ description (UG): White, round, flat tablets embossed "UG" above and "KPI " below the break-line on one side and double cross line on the other.		Upper:	Flat Beveled, Embossed UG over B/L and KPI under B/L
		Lower:	Flat, beveled with double cross line
Cleaned by:		Date:	
Machine setting done by:		Date:	
Machine speed:			
Checked by (section supervisor) :		Date:	
Previous Product:		Batch No.:	

3. Tablet Weight Range Limits

Individual tablet wt. (Nominal):			674.6mg
% of Nominal weight		Tablet wt. (mg)	Symbol
Upper range	>105%	>708.3	> +T2 Alert
	105%	708.3	= +T2 Alert
	103%	694.8	= +T1 Good
Norm	100%	674.6	Very Good
Lower range	97%	654.4	= -T1 Good
	95%	640.9	= -T2 Alert
	<95%	<640.9	< -T2 Alert

4. COMPRESSION

Started	Date:	Time:	Completed	Date:	Time:
---------	-------	-------	-----------	-------	-------

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 35 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
TABLET COMPRESSION				

5. DETERMINE PERCENTAGE YIELD

Date:	Shift:	Operator/ Grp. Leader	Weight of tablets compressed
Total weight of compressed tablets			Kg

YIELD RECONCILIATION

Steps	Descriptions	Qty. (Kg)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size with respect to Blended Granules			
B.	Actual Quantities of Tablets Compressed			
C.	Sample Quantities (In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
Percentage Yield =F x100/A =.....%				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):
CAUSE OF VARIATION: YES/NO

REMARKS/BREAKDOWNS.....

.....
.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 36 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
TABLET COMPRESSION				

6. INITIAL COMPRESSION MACHINE SETTING INSTRUCTIONS:

NOTE: Collect the tablets of first three rounds and discard as non-recoverable residues.

1. Run feeder for one minute.
2. Set the weight of tablets as per limits given above by turning weight dossier anticlockwise for increasing weight and clockwise for decreasing weight & discard the tablets.
3. Check the weights of 20 tablets in group. Set weights by turning the weight dossier anticlockwise for increasing the weight and clockwise for decreasing the weight. Check total weights of twenty tablets three times & record. Calculate Average weight. Check the individual weights of all punches and record.
4. Initially check thickness of all stations tablets each side and during process.
5. Initially check hardness of all stations tablets each side and during process.
6. Check friability, disintegration time and individual weight variation of tablets.
7. Collect fresh samples after complete setting of parameters and visually inspect for any defects.
8. If it is found satisfactory then proceed for further process and verify all in process parameters and record as per next table.

STANDARD PARAMETERS		
PARAMETERS	FREQUENCIES	QUANTITIES
Appearance	Initial	All Station Tabs.
	Every 60 Minutes	20 Tablets/Side
Marking/ Embossing	Initial & Every 60 Minutes	All Station Tab.
Group Weight	Initial & Every 60 Minutes	20 Tablets/Side
Individual weight	Initial	All Station Tabs.
	Every 60 Minutes	20 Tablets/Side
Disintegration Time	Initial & Every 60 Minutes	6 Tablets/Side
Uniformity	Every 60 Minutes	-----
Hardness	Initial	All Station Tabs.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 37 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
TABLET COMPRESSION				

	Every 60 Minutes	All Station Tabs.
Thickness	Initial	All Station Tabs.
	Every 60 Minutes	All Station Tabs.
Friability	Every 60 Minutes	20 Tabs/side
M/c RPM	-----	-----

SPECIFICATIONS:

S.N.	TESTS	SPECIFICATION
1.	Physical Appearance:	White round flat tablets embossed 'F' on one side and crossed break lines on the other.
2.	Diameter Limits (mm)	12.5 ± 0.1
3.	Thickness Limits (mm)	4.48 ± 0.3
4.	Hardness Limits (Kg/cm ²)	4-12
5.	Tablet Standard Weight (mg)	674.6 ± 2%
6.	Weight Uniformity Limits	± 5% of Average Weight. NMT 2 tablets should be out of 5% range and none of the tablets should be out of 10% range.
7.	Friability Limits	NMT 1.0% If obviously cleaved, cracked, chipped or broken tablets are present in the sample after tumbling, the sample fails the test.
8.	Disintegration time (Min)	NMT 15.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 38 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

7. INITIAL DIES AND PUNCHES CHECKING RECORD

STA. No.	UPPER PUNCH	LOWER PUNCH	DIES SUITABLE	STA. No.	UPPER PUNCH	LOWER PUNCH	DIES SUITABLE	STA. No.	UPPER PUNCH	LOWER PUNCH	DIES SUITABLE
01				21				41			
02				22				42			
03				23				43			
04				24				44			
05				25				45			
06				26				46			
07				27				47			
08				28				48			
09				29				49			
10				30				50			
11				31				51			
12				32				52			
13				33				53			
14				34				54			
15				35				55			
16				36				56			
17				37				57			
18				38				58			
19				39				59			
20				40				60			

PUNCH SET BY (Operator Name)

CHECKED BY (Production Supervisor)

VERIFIED BY(QA).....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 39 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

8. INITIAL INDIVIDUAL WEIGHT RECORDS

Station. Nos.	Individual weights	Station. Nos.	Individual weights	Station. Nos.	Individual weights	Station. Nos.	Individual weights	Station. Nos.	Individual weights
LHS									
01		12		23		34		45	
02		13		24		35		46	
03		14		25		36		47	
04		15		26		37		48	
05		16		27		38		49	
06		17		28		39		50	
07		18		29		40		51	
08		19		30		41		52	
09		20		31		42		53	
10		21		32		43		54	
11		22		33		44		55	
RHS									
01		12		23		34		45	
02		13		24		35		46	
03		14		25		36		47	
04		15		26		37		48	
05		16		27		38		49	
06		17		28		39		50	
07		18		29		40		51	
08		19		30		41		52	
09		20		31		42		53	
10		21		32		43		54	
11		22		33		44		55	
Done by Operator (Sign & Date)			Checked by Supervisor (Sign & Date)			Verified by QA (Sign & Date)			

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 40 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

9. TABLET COMPRESSION IN-PROCESS QUALITY CONTROL REPORT									
M/C No.:				Location				Speed	
(I) Tablet Weight In-Process Control									
Time HH:MM									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance									
Individual Weights	1								
	2								
	3								
	4								
	5								
	6								
	7								
	8								
	9								
	10								
	11								
	12								
	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
Total wt. (mg)									
Wt. after Friability									
Avg wt. (mg)									
Range									
Max. wt. (mg)									
Min. Wt. (mg)									
(II) Max/min wt. plot									
Action >+T2	>708.3mg								
Alert+T2	708.3mg								
Good +T1	694.8mg								
V. Good Norm	674.6mg								
Good -T1	654.4mg								
Alert -T2	640.9mg								
Action <-T2	<640.9mg								
Friability	NMT 1%								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 41 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

(III) Other In-Process Control Tests									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Disintegration	NMT 15min								
Hardness	4-12 (kg/cm ²)								
Diameter (mm)	12.50 ± 0.1								
Thickness (mm)	4.48 ± 0.3								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 42 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

TABLET COMPRESSION IN-PROCESS QUALITY CONTROL REPORT									
M/C No.:				Location				Speed	
(I) Tablet Weight In-Process Control									
Time HH:MM									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance									
Individual Weights	1								
	2								
	3								
	4								
	5								
	6								
	7								
	8								
	9								
	10								
	11								
	12								
	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
Total wt. (mg)									
Wt. after Friability									
Avg wt. (mg)									
Range									
Max. wt. (mg)									
Min. Wt. (mg)									
(II) Max/min wt. plot									
Action >+T2	>708.3mg								
Alert+T2	708.3mg								
Good +T1	694.8mg								
V. Good Norm	674.6mg								
Good -T1	654.4mg								
Alert -T2	640.9mg								
Action <-T2	<640.9mg								
Friability	NMT 1%								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 43 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

(III) Other In-Process Control Tests									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Disintegration	NMT 15min								
Hardness	4-12 (kg/cm ²)								
Diameter (mm)	12.50 ± 0.1								
Thickness (mm)	4.48 ± 0.3								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 44 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

TABLET COMPRESSION IN-PROCESS QUALITY CONTROL REPORT									
M/C No.:				Location				Speed	
(I) Tablet Weight In-Process Control									
Time HH:MM									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance									
Individual Weights	1								
	2								
	3								
	4								
	5								
	6								
	7								
	8								
	9								
	10								
	11								
	12								
	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
Total wt. (mg)									
Wt. after Friability									
Avg wt. (mg)									
Range									
Max. wt. (mg)									
Min. Wt. (mg)									
(II) Max/min wt. plot									
Action >+T2	>708.3mg								
Alert+T2	708.3mg								
Good +T1	694.8mg								
V. Good Norm	674.6mg								
Good -T1	654.4mg								
Alert -T2	640.9mg								
Action <-T2	<640.9mg								
Friability	NMT 1%								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 45 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

(III) Other In-Process Control Tests									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Disintegration	NMT 15min								
Hardness	4-12 (kg/cm ²)								
Diameter (mm)	12.50 ± 0.1								
Thickness (mm)	4.48 ± 0.3								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 46 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

TABLET COMPRESSION IN-PROCESS QUALITY CONTROL REPORT									
M/C No.:				Location				Speed	
(I) Tablet Weight In-Process Control									
Time HH:MM									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance									
Individual Weights	1								
	2								
	3								
	4								
	5								
	6								
	7								
	8								
	9								
	10								
	11								
	12								
	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
Total wt. (mg)									
Wt. after Friability									
Avg wt. (mg)									
Range									
Max. wt. (mg)									
Min. Wt. (mg)									
(II) Max/min wt. plot									
Action >+T2	>708.3mg								
Alert+T2	708.3mg								
Good +T1	694.8mg								
V. Good Norm	674.6mg								
Good -T1	654.4mg								
Alert -T2	640.9mg								
Action <-T2	<640.9mg								
Friability	NMT 1%								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 47 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

(III) Other In-Process Control Tests									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Disintegration	NMT 15min								
Hardness	4-12 (kg/cm ²)								
Diameter (mm)	12.50 ± 0.1								
Thickness (mm)	4.48 ± 0.3								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 48 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

TABLET COMPRESSION IN-PROCESS QUALITY CONTROL REPORT									
M/C No.:				Location				Speed	
(I) Tablet Weight In-Process Control									
Time HH:MM									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Physical Appearance									
Individual Weights	1								
	2								
	3								
	4								
	5								
	6								
	7								
	8								
	9								
	10								
	11								
	12								
	13								
	14								
	15								
	16								
	17								
	18								
	19								
	20								
Total wt. (mg)									
Wt. after Friability									
Avg wt. (mg)									
Range									
Max. wt. (mg)									
Min. Wt. (mg)									
(II) Max/min wt. plot									
Action >+T2	>708.3mg								
Alert+T2	708.3mg								
Good +T1	694.8mg								
V. Good Norm	674.6mg								
Good -T1	654.4mg								
Alert -T2	640.9mg								
Action <-T2	<640.9mg								
Friability	NMT 1%								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 49 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET COMPRESSION					

(III) Other In-Process Control Tests									
Done by:		Operator:		QA:		Operator:		QA:	
Tablet discharge		Left	Right	Left	Right	Left	Right	Left	Right
Disintegration	NMT 15min								
Hardness	4-12 (kg/cm ²)								
Diameter (mm)	12.50 ± 0.1								
Thickness (mm)	4.48 ± 0.3								

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 50 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
TABLET COMPRESSION				

10. BULK TRANSFER / RECONCILIATION
From: Tablet compression Section
To: Semi-Finished product quarantine store (SMFPQ)

Please receive the following tablets:

Date	Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)	Delivered by Operator	Received by SMFPQ Supervisor
	1					
	2					
	3					
	4					
	5					
	6					
	7					
	8					
	9					
	10					
Total						
Total number of tablets:				Average weight of tablets:		

Section supervisor (sign & date) _____

Remarks

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 51 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
TABLET INSPECTION & SORTING					

LINE CLEARANCE FORM					
Previous product:		Batch No:			
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign &	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. Area cleaning as per SOP.					
2. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.					
3. Gowning procedure is followed.					
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.					
5. Check the proper cleanliness of the Air Supply and Return grills and Air Conditioning System.					
6. Ensure that the Vibrosifter sieves, Metal Detector, Deduster, Scoops etc. are cleaned.					
7. Check and ensure that Balance Calibration & Verification records are updated.					
8. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (40 to 65%) _____					
Remarks:					
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign &	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
1. All material from current product have been removed.					
2. All facilities have been cleaned and labeled.					
3. The area has been cleaned and labeled.					
4. All paperwork for the current batch has been completed.					
5. Reconciliation of current batch has been done.					
6. Machine(s) has been checked.					
Remarks:					

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 52 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
TABLET INSPECTION & SORTING				

INSPECTION AND SORTING PROCEDURE

1. Ensure that the inspection area and equipment are clean.
2. Ensure that all personnel are properly dressed in protective garments, gloves and masks.
3. Carry out a line clearance exercise as per SOP for line clearance and complete the "beginning activities" section of the record.
4. Bring the tablets inside the sorting room.
5. Visually inspect the tablets for defects such as chipping, spots etc.
6. Collect all rejected tablets in properly labeled designated containers.
7. Seal inspected tablets in polythene lined drums and properly label the drums.
8. Complete the "Ending activities" section of the line clearance form.
9. Transfer inspected and sorted tablets to Packing area/Semi-finished product quarantine (SMFPQ) store.

VISUAL INSPECTION LOG

ACTIVITY	FROM	TO

RECONCILIATION AFTER INSPECTION AND SORTING:

1	Weight of tablets received for inspection and sorting _____ Kg
2	Weight of tablets received after inspection and sorting _____ Kg
3	Average weight of tablets _____ mg/tablet
4	Recyclables after inspection and sorting _____ Kg (recyclables are recorded in the BMR reconciliation sheet)
5	Total rejects of tablets after inspection and sorting _____ Kg (rejects are recorded in the BMR reconciliation sheet).
6	Total loss after inspection and sorting _____ Kg _____ %

Done by **group leader (sign & date):**

Checked by **SMFPQ supervisor (sign & date):**

Verified by **QA (sign & date):**

7	Names of person(s) involved in carrying out Inspection and Sorting		
	(a)	(b)	(c)
	(d)	(e)	(f)
	(g)	(h)	(i)

NOTE: Transfer inspected and sorted tablets to Semi Finished Product Quarantine (SMFPQ) Store.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 53 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
TABLET INSPECTION & SORTING				

QA IN-PROCESS QUALITY CONTROL REPORT FOR INSPECTION & SORTING
--

Product compressed from Compression Machine Number _____

Record the net content of the drums to be sorted, weigh and record the quantity of defects sorted from the respective drums in the table below.

Drum number	Net weight (kg)	Quantity of defects (kg)	Nature of defects
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			
Total			
Percentage defects = (Total quantity of defects / total net weight) x 100% = _____ %			

Done by (QA officer, sign & date): _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 54 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	

PACKAGING MATERIALS REQUISITION SHEET

Requested by Production Supervisor		Name:		Signature:		Date:					
Approved by Production Manager		Name:		Signature:		Date:					
Pack Size: (/ mark as appropriate)				10x10 Blisters	100 Blisters X 10 Tablets	500's	1000's				
S/N	Item Code	Item Description	AR. NO	Total Batch Quantity	Units	Issued by (PM Stores In-charge)		Received by (Packaging Supervisor)		Verified by (QA)	
						Sign	Date	Sign	Date	Sign	Date
FOR BLISTER PACK											
1	2KSD2TB	Kamsidar Tablets Aluminium Foil 203 x 0.25 mm (D/T)			Kg						
2	1KSD1TB	Kamsidar Tablets Aluminium Foil 205 x 0.25 mm (D/T)			Kg						
3	4KSD4TB	Kamsidar Tablets Aluminium Foil 205 x 0.25 mm (D/T)-UG									
4	PF 206M	PVC Clear Film 206 x 0.25mm (D/T)			Kg						
5	PF 210M	PVC Clear Film 210 x 0.25mm (D/T)			Kg						
6	1KSD1TB	Kamsidar Tablets 10 x 10 blister cartons			PCS						
7	STB1XT	350 x 10 x 10's Tablets blister shipping cartons			PCS						
8	BPTPE2	2" BOPP Printed Packing Tape			PCS						
9	PILKSD	Kamsidar patient information leaflet			PCS						
FOR BULK PACK											
1	3KSD4TJ	Kamsidar Tablets 1000's Labels			PCS						
2	3KSD4TJ	Kamsidar Tablets 1000's Labels (UG)			PCS						
3	PJ1KCC	1000cc Plastic Jars With Caps			PCS						
4	LDB810	8"x10" HDPE Plain, Clear Bags			PCS						
5	SJJ4XM	20x1000cc Shipping Cartons			PCS						
6	BPTPE2	2" BOPP Printed Packing Tape			PCS						
7	PILKSD	Kamsidar patient information leaflet			PCS						
8	3KSD4TJ	Kamsidar Tablets 1000's Labels			PCS						

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 55 of 72	
BATCH SIZE		600,000 Tablets	Mfg. Date :		
BATCH No.			Exp. Date :		
MFR No.	KPI/MFR/022/00	Date:	Shift :		
PACKAGING					

LINE CLEARANCE FORM				
Previous product:		Batch No:		
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date) Status
1. Area cleaning as per SOP.				
2. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.				
3. Gowning procedure is followed.				
4. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.				
5. Check the proper cleanliness of the Air Supply and Return grills and Air Conditioning System.				
6. Ensure that the Vibrosifter sieves, Metal Detector, Deduster, Scoops etc. are cleaned.				
7. Remove all the cartons, labels, strips, any material left over on the line and on the lower space of the conveyor belts.				
8. Put an authorized overprinted specimen of strip/blister/cartons/Label/Tube in the BMR for ready reference.				
9. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (40 to 65%) _____				
Remarks:				
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by QA (Sign & Date) Status
1. All material from current product have been removed.				
2. All facilities have been cleaned and labeled.				
3. The area has been cleaned and labeled.				
4. All paperwork for the current batch has been completed.				
5. Reconciliation of current batch has been done.				
6. Stereos of previous product are submitted to production Supervisor.				
Remarks:				

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 56 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
PACKAGING					

1. PACKING PROCESS & INSTRUCTIONS

BULK PACKING

1. Ensure that the packing room, packing equipment and packing table(s) are clean.
2. Carry out a Packing line clearance exercise as per SOP for line clearance and complete the "beginning activities" section of the record.
3. Bring in the appropriate packaging material.
4. Check the labels for correctness and legibility of coding.
5. Bring in the tablets to be packed.
6. Count/Weigh the tablets.
7. Put the tablets in poly bags and seal.
8. Place the sealed poly bag in the container and firmly place the cap.
9. Re-check the batch no., mfg. date and exp. date on the labels.
10. Carefully label the container.
11. Pack the right number of containers in the appropriate corrugated shipping cartons.
12. Seal the shipping cartons with packing tape.
13. Label the shipping cartons appropriately.
14. Complete the packing records.
15. Complete the "Ending activities" section of the packing line clearance form.
16. Transfer packed goods to Finished Goods Quarantine Store and wait for QC report.

BLISTER PACKING

1. Ensure that the packing room, packing equipment and packing table(s) are clean.
2. Carry out a Packing line clearance exercise as per SOP for line clearance and complete the "beginning activities" section of the record.
3. Bring in the appropriate packaging material.
4. Set and operate the blister packing machine as per SOP for operating blister packing machine.
5. Check the first code for correctness and legibility of coding.
6. Bring in the tablets to be blistered and complete the operation and transfer all those blisters to secondary packing area.
7. Pack 10 blisters per unit carton/ 100 blisters per unit dispenser.
8. Pack the right number of unit cartons in the corrugated shipping cartons.
9. Seal the shipping carton with packing tape.
10. Label the shipping carton appropriately.
11. Complete the packing records.
12. Complete the "Ending activities" section of the packing line clearance form.
13. Transfer packed goods to Finished Goods Quarantine Store and wait for QC report.

Reference Documents:

1. SOP for status labelling.
2. SOP for line clearance.
3. SOP for Bulk packing.
4. SOP for Blister packing

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 57 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
PACKAGING					

2. TABLET BLISTERING

1. BLISTER MACHINE SETTING

Blister Machine number		Blister machine Make	
PVC Clear Film Size		Aluminium Foil Size	
Set Sealing Temperature		Observed Sealing Temperature	
Set Forming Temperature		Observed Forming Temperature	
Machine Speed		Blister Pack Size	
Machine Setting Done By (Sign & Date):			
Checked by Supervisor (Sign & Date):			
Verified by QA (Sign & Date):			

2. BLISTER PACKING

Started	Date:	Time:	Completed	Date:	Time:
---------	-------	-------	-----------	-------	-------

3. DETERMINE PERCENTAGE YIELD

Date:	Shift:	Operator/ Grp. Leader	Number of Blisters
Total Number of Blisters			

YIELD RECONCILIATION

Steps	Descriptions	Qty. (units)	Done by SPV (Sign & Date)	Verified by QA (Sign & Date)
A.	Theoretical Batch Size (Blisters) with respect to tablets received from sorting			
B.	Actual Quantities of Blisters			
C.	Sample Quantities (In-Process checks)			
D.	Validation Samples			
E.	Rejects (if any)			
F.	Total Actual Yield (B+C+D)			
G.	Unaccountable Losses [A-(E+F)]			
Percentage Yield = $F \times 100 / A = \dots\dots\dots\%$				
Permissible Yield Limit: (98% to 102%)				

(If percentage yield is within 98% to 102% no justification is required, otherwise explain variation in space below):

CAUSE OF VARIATION: YES/NO

REMARKS/BREAKDOWNS.....

.....

.....

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 58 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
PACKAGING					

3. PACKING CODING CONTROL

1	Batch No:	Coding approved by Packing Supervisor (sign & date)	
2	Mfg. Date:		
3	Exp.Date:		
4	Item(s) to be coded: Labels, Unit Cartons, Shipper labels, Blisters, Others (specify); mark as appropriate		
5	Coding Stamp set by (operator): sign & date:		
	Coding Machine Number:		
6	Attach first coded samples and patient information leaflet (PIL) here:		
7	Coding approved by: QA Officer (Sign & Date):		
8	Names of person(s) applying labels and packing - BULK & BLISTER (mark as appropriate).		
	(a)	(b)	(c)
	(d)	(e)	(f)
	(g)	(h)	(i)
9	Coding Reconciliation:		
	Name of item:		
	No. of item issued:		
	No. of item coded:		
	No. of item uncoded:		
	No. of items damaged / rejected during coding:		
	No. of items used:		
	No. of excess items, if any:		
	No. of item to be destroyed (with QA consent):		
	Reconciliation approved by: QA (Sign & Date):		

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 59 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
PACKAGING					

4. BULK TRANSFER AFTER INSPECTION & SORTING

From: Semi finished product quarantine store (SMFPQ)

To: Packing Department (Bulk / blisters/strips)

Please receive the following tablets for packing:

Date	Drum No.	Gross Weight (Kg)	Tare Weight (Kg)	Net Weight (Kg)	Delivered by SMFPQ supervisor	Received by Operator/ Packing supervisor
	1					
	2					
	3					
	4					
	5					
	6					
	7					
	8					
	9					
	10					
Total						
Total number of tablets:				Average weight of tablets:		

Quantity packed (bulk /blisters): -----

(A) REJECTS AFTER PACKING _____ KG/G

(B) RE-WORKABLES AFTER PACKING _____ KG/G

Section supervisor (sign & date) _____

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 62 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
PACKAGING				

[illegible]

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 63 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SECONDARY PACKAGING					

LINE CLEARANCE FORM					
Previous product:		Batch No:			
BEGINNING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements.		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
10. Area cleaning as per SOP.					
11. Ensure that all the previous product containers, material, and labels are removed from the manufacturing area.					
12. Gowning procedure is followed.					
13. Ensure that status board is displayed mentioning Product Name, Batch No, B. Size, Mfg. Date, Exp.Date & Status with sign & date.					
14. Check the proper cleanliness of the Air Supply and Return grills and Air Conditioning System.					
15. Ensure that the Vibrosifter sieves, Metal Detector, Deduster, Scoops etc. are cleaned.					
16. Remove all the cartons, labels, strips, any material left over on the line and on the lower space of the conveyor belts.					
17. Put an authorized overprinted specimen of strip/blister/cartons/Label/Tube in the BMR for ready reference.					
18. Ensure Environmental conditions are met as per SOP. Temperature NMT 28 °C _____ Relative Humidity (40 to 65%) _____					
Remarks:					
ENDING ACTIVITIES: The following must be checked to ensure adherence to GMP requirements		Done by Operator (Sign & Date)	Checked by Supervisor (Sign & Date)	Verified by	
				QA (Sign & Date)	Status
7. All material from current product have been removed.					
8. All facilities have been cleaned and labeled.					
9. The area has been cleaned and labeled.					
10. All paperwork for the current batch has been completed.					
11. Reconciliation of current batch has been done.					
12. Stereos of previous product are submitted to production Supervisor.					
Remarks:					

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 64 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.		KPI/MFR/022/00		Date:	
				Shift :	
SECONDARY PACKAGING					

INSPECTION AND SORTING OF BLISTERS AND STRIPS			
INSPECTION AND SORTING PROCEDURE			
- All personnel involved in Inspection, Sorting and Packing shall be qualified annually with eye vision checks (Ocular Test), colour-blindness checks, and ability to distinguish known product rejects from acceptable products. - The medical examination of all Inspectors/Sorters and Packers should include the report of proper vision with or without glasses and Color blindness. - The production Inspectors/Sorters and Packers who feel uneasy or complain of headache are unfit for handling Inspection, Sorting and Packing of Blisters/Strips. - Inspectors/Sorters who are employed for checking strips and blisters or overprinted materials shall be replaced with other trained checkers alternately after every four hours or a break of at least 30 minutes shall be taken to counter eye fatigue and complacency. - Before the start of any batch, line clearance should be done as per SOP for line clearance. - Bring the strips/blisters to the secondary packing room. - Visually inspect/Check the strips / blisters / coded materials for all the possible defects. Some examples of types of rejection on packing line are Coding rejection, Due to text matter/hologram/smudging of ink over printed foil, Rejected strips due to cutting, Empty pocket rejection, Pocket puncture rejection, Sealing rejection, Overheated strips/blisters, Chipped / capped / broken tablets, Improper forming of blisters, Improper appearance / absence of hologram in the final blister strip among others. - After the Inspection and Sorting operation, the Supervisor shall intimate to the QA Officer online to verify the integrity of the sorted blisters/strips and approve packing to start. - Remove all defective blisters/strips, place them in a polybag or drum and label as recyclables. These should be de-blistered/de-stripped once sorting and inspection of blisters is done.			
Inspection and Sorting Line Clearance Done by (sign & date)			
Checked by Production Supervisor (sign & date):			
Verified by QA (sign & date):			
VISUAL INSPECTION LOG			
ACTIVITY	FROM	TO	
10	Names of person(s) involved in carrying out Inspection and Sorting		
	(a)	(b)	(c)
	(d)	(e)	(f)
	(g)	(h)	(i)

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 65 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SECONDARY PACKAGING					

SECONDARY PACKING PROCEDURE			
1. Check the specimen attached to the BMR of overprinted carton/Monocartons, leaflet, and shipper label if it has been checked by the line supervisor and approved by the QA Officer online. 2. Transfer previously approved overprinted packing material from the printing room to the secondary packing room. 3. Check if displayed details on the carton/Monocartons/dispensers of the batch number, manufacturing and expiry date are as per BMR/BPR. Check carton for overprinting intactness. 4. Check if the information on leaflets/inserts is clear, legible and specific to the product to be packed. 5. Open the carton and insert leaflet with number of strips/blisters as specified in the BMR. Close the carton and arrange it into the shipper. 6. Once the shipper is filled with the number of cartons as specified in the BMR, insert a fully filled packing slip with name of the product & batch number, the name of the packing personnel and this should be checked and signed by the group leader. 7. Group leaders should check the shipper to ensure that it has the correct quantity of cartons/Monocartons/dispensers. If this is okay, group leader shall sign the packing slip and intimate to the Supervisor and QA Officer online to confirm the correctness of the quantity in the shipper and other packing requirements. The Supervisor and the QA officer shall then sign on the packing slip and authorize the group leader to seal the shipper with BOPP tape and attach the shipper label. Exact quantity of loose shipper is written on shipper label and sealed. 8. Once the entire batch has been completed, transfer to finished product quarantine store. 9. Any leftover printed cartons/Monocartons/dispensers/labels should be destroyed in presence of QA Officer online and the number recorded in the BMR. 10. Transfer leftover packaging material (that is not printed on) to Packaging Material Stores with a return note.			
Secondary Packing Line Clearance Done by (sign & date)			
Checked by Production Supervisor (sign & date):			
Verified by QA (sign & date):			
RECONCILIATION AFTER PACKING			
Product Name:		Batch Number:	
Batch Size (Cartons/Jars):		Pack Size:	
No. of intact shippers packed		No. of cartons/jars in intact Shippers	
Loose Quantity Packed:		Total Quantity Packed (cartons/jars)	
Done by group leader (sign & date):			
Checked by SMFPQ supervisor (sign & date):			
Verified by QA (sign & date):			
11	Names of person(s) applying labels and packing - BULK & BLISTER (mark as appropriate).		
	(a)	(b)	(c)
	(d)	(e)	(f)
	(g)	(h)	(i)

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 66 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
SECONDARY PACKAGING				

SHIPPER WEIGHT VERIFICATION

- After the packing operation is completed for the batch, **Weigh the packed shippers one by one** and record the weight of packed shipper in BPR and also on the Shipper Label and recheck the shipper if any discrepancy is found.

Average Weight of Empty Carton:					Average Weight of Packed Unit Carton:						
Average Weight of Empty Shipper:					Average Weight of Packed FULL Shipper:						
Acceptance Limit: _____ % of average weight of packed Shipper					Range: _____ to _____						
SHIPPER No.	SHIPPER WEIGHT (Kg)	SHIPPER No.	SHIPPER WEIGHT (Kg)	SHIPPER No.	SHIPPER WEIGHT (Kg)	SHIPPER No.	SHIPPER WEIGHT (Kg)	SHIPPER No.	SHIPPER WEIGHT (Kg)	SHIPPER No.	SHIPPER WEIGHT (Kg)
01		11		21		31		41		51	
02		12		22		32		42		52	
03		13		23		33		43		53	
04		14		24		34		44		54	
05		15		25		35		45		55	
06		16		26		36		46		56	
07		17		27		37		47		57	
08		18		28		38		48		58	
09		19		29		39		49		59	
10		20		30		40		50		60	
Remarks: All shippers are within Specifications YES/NO _____											
Weighing Done by (Sign & Date):											
Checked by Packing Supervisor (Sign & Date):											
Approved by QA (Sign & Date):											

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 67 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
SECONDARY PACKAGING					

[illegible]

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR		Page 69 of 72
BATCH SIZE		600,000 Tablets	Mfg. Date :	
BATCH No.			Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:	Shift :	
SECONDARY PACKAGING				

[illegible]

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 70 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Shift :	
FINISHED PRODUCT RECONCILIATION SHEET					

Quantity Packed	Retention Sample Quantity	Actual Quantity Released	Quantity Delivered to FGS	Delivery Note Number	Delivery Date

FINAL BATCH RECONCILIATION

ITEM	Quantity (Gm/Kg).	Quantity* (No. of Units)	Percent =[QTYx100/(A)]
(A) Theoretical Yield =(Batch Size):			
(B) Yield - 1st Crop (Quantity Packed):			
(C) In-Process/QC Samples:			
(D) Re-workables =[(I)+(II)]:			
(I) Powder/Granules*			
(II) Tablets			
(E) Rejects =[(I)+(II)]:			
(I) Powder/Granules*			
(II) Tablets			
(F) Total Yield =(B+C+D):			
(H) Shortage =[A-(E+F)]:			

Done by Supervisor (Sign & Date):	Verified by QAO (Stamp, Sign & Date):
--	--

Reconciliation Approved By **Production Manager (Sign & Date)**: _____

*Quantity (No. of units) is derived by converting this quantity using average compression weight.

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.	
PRODUCT		KAMSIDAR		Page 71 of 72	
BATCH SIZE		600,000 Tablets		Mfg. Date :	
BATCH No.				Exp. Date :	
MFR No.	KPI/MFR/022/00	Date:		Actual Quantity :	
BATCH RELEASE SHEET					

The following Batch production records have been properly filled out, checked, dated and signed by authorized person(s). The tests are also within limits. The batch is hereby approved for release.		
Reviewed and approved by QA Head/Designee (Sign & Date):		
Records Reviewed (Mark as appropriate)	QA Head/Designee	Company /Prod. Pharmacist
1. Raw Material Requisition/Dispensing Sheet		
2. Mixing / Granulation Records		
3. Blending / Lubrication Records		
4. Compression/ In-process Control Records		
5. Process Line Clearance Records		
6. Records for Semi-finished Products (blended material, compressed tablets etc.)		
7. Inspection and Sorting Records		
8. Packing Line Clearance Records		
9. Packing Material Authorization Note		
10. Packing Coding Control Records		
11. Product Reconciliation Sheet		
12. IPQA Reports		
13. QC Analytical reports		
13.1. Finished Product Analytical Report		
13.2. Finished Product Analytical Data Sheet (Observations, Calculations, Printouts, References).		
I hereby certify that all the manufacturing stages of this batch of finished product have been carried out in full compliance with the Uganda National Drug Authority (NDA) GMP requirements and with the requirements of the product registration. I, therefore, AUTHORISE release of the batch.		
Company Pharmacist/Pharmacist In-Charge:		
_____	_____	_____
Name	Signature	Date

		Batch Manufacturing Record and Procedure		KAMPALA PHARMACEUTICAL INDUSTRIES (1996) LTD.
PRODUCT		KAMSIDAR	Page 72 of 72	
BATCH SIZE		600,000 Tablets		
MFR No.	KPI/MFR/022/00			
REVISION HISTORY				

REVISION NO.	CHANGES INCORPORATED	REASON OF CHANGE	EFFECTIVE DATE
01	- KPI new logo to be incorporated - The provision for BMR preparation, Review, Approval and Authorization have been introduced on the Master Formula Page(s). - An Index has been introduced. - The Abbreviations used in the BMR have been defined in full. - List Of Equipment & Instruments used in product manufactured have been summarized. - Line Clearance for Dispensing previously missing has been introduced. - Documentation of First Drying, Second Drying, and Final Drying have been split and captured independently, and Milling activities including Sieves used plus LOD data are also being documented at granulation. Identity of the FBD used to dry a particular Lot has been captured. - In-process Quality Control and Assurance (IPQC/IPQA) reports for both Production and QA at all stages (as appropriate) have been incorporated. - Initial Compression Machine Setting instructions, Initial Dies and Punches Checking Record, Initial Individual Weight Records as per number of Punches/Stations of the Compression machine have been introduced/Captured.	To improve on clarity of the BMR. To match up with the current practices.	
02	- To eliminate duplication in the BMR header, introduce page breaks for different process steps and include MFR numbers in the BMR. - To include the Line Clearance for Inspection and Sorting of Blisters/Strips and the Secondary Packing Procedure.	To improve coherence of the BMR To ensure maximum elimination of blister/Strip defects before packing.	
03	- Batch Release sheet revised to ensure that batch release is Authorized by the Company Pharmacist/Pharmacist In-Charge. - Shipper weight verification record introduced. - Packaging Material Requisition Sheet incorporated in the BMR.	To comply with regulatory requirements for batch release. To eliminate the risk of having less or more product quantity in a shipper	